



Quantum IV

Body Composition Analyzer



Measuring and Managing Healthy Living

rjlsystems.com

Quantum IV is our flagship product and represents a culmination of more than 30 years of technology and continuous improvement. It is lightweight, easy to use, and fully portable. The instrument case has been designed to withstand almost any environment, ensuring the accuracy and reliability of your Quantum IV device.

Key Features

The Quantum IV features a bright green LED display that is easy to read and provides all the pertinent information on one screen in your choice of English, Spanish, Italian, or German. Not only does this simplify the testing procedure, but it allows the entire test to be conducted very quickly.

Affordable

This professional Quantum IV BIA Analyzer is affordable, effective and accurate.

Easy To Read

Bright green LED display shows resistance, reactance, impedance.

Software

RJL Systems BC 4.0 Software is Windows and Mac Compatible.

Portable


The Quantum IV weighs just 860 grams and has a comfortable ergonomic design.

Powerful

Rechargeable lithium-ion batteries run for hours, making your device portable and ready to go.

Update Ready

As RJL develops new technologies, the Quantum IV and corresponding software will include new features and firmware.



Quantum IV has the ability to run for hours continuously, before a simple recharge of its internal lithium ion batteries, so you can be assured that your Quantum IV is always ready to go. Adding to the legendary quality of the RJL Quantum BIA, this robust feature set also includes an ergonomically designed case, ensuring years of reliable service even in the harshest environments.

Included:	Quantum IV	Quantum IV+
Quantum IV Body Composition Analyzer (\$3,090 ⁰⁰)	✓	✓
BC 4.0 Software (Mac and PC) (\$500 ⁰⁰)	✓	✓
Quantum IV Warranty (\$99/year)	✓	✓✓✓
Software upgrade protection and telephone support (\$200 ⁰⁰ /year)	✓	✓✓✓
Subject cables (\$60 ⁰⁰ /per set)	✓	✓✓✓
400 electrodes (\$100 ⁰⁰)	✓	✓✓✓✓✓ ✓✓✓✓✓*
Technician Webinar (Complimentary)	✓	✓
Lifestyle Practitioner Webinar (Professionally Recorded) (\$99 ⁰⁰)		✓
Advanced Practitioner Webinar (Live) (\$149 ⁰⁰ /\$169 ⁰⁰)		✓
Value	\$4,049⁰⁰	\$5,315⁰⁰
Price	\$3,090⁰⁰	\$3,990⁰⁰
	YOU SAVE \$959⁰⁰	\$1,325⁰⁰

* E-Coupons to ensure freshness (a \$900 value!)